

ML COURSE

NTI

**Intelligent Multi-Condition
Medical Prediction System**

PROJECT OVERVIEW

Our system handles three different medical prediction tasks:

Medical Condition	Data Type	Prediction Task	Output
 ECG	ECG Signals	Multiclass Classification	Heartbeat Class
 Heart Failure	Clinical Data	Binary Classification	Death / No Death
 CKD	Clinical Data	Binary Classification	CKD / No CKD

ONE AI SYSTEM → THREE MEDICAL PREDICTION TASKS

Problem Statement

◀ Medical data can be:

- Large and complex
- Difficult to analyze manually
- Time-consuming to interpret
- Rich in hidden patterns

◀ Our Solution

Use AI to learn patterns from medical data and provide predictions.



PROJECT OBJECTIVES



- Analyze three medical datasets
- Perform EDA and data cleaning
- Apply appropriate preprocessing techniques
- Train multiple ML and DL models



- Evaluate and compare model performance
- Select the best model for each task
- Deploy the models into one web application

DATASETS



Dataset	Data Type	Task	Target
ECG	Signal	Multiclass Classification	5 Classes
Heart Failure	Tabular	Binary Classification	DEATH_EVENT
CKD	Tabular	Binary Classification	CKD / No CKD

- ◀ **Signal = data representing a changing signal over time**
- ◀ **Tabular = data organized into rows and columns**

Project Workflow:



MEDICAL DATA INPUT (Telemetry & Tabular Records)



ECG Pipeline

Preprocessing (Scaling & R-Peak Filtering)

Model: Gradient Boosting / RF

Evaluation & Waveform Plotting

Heart Failure

Preprocessing (Selective Scaling)

Model: Survival Classifier

Evaluation & Risk Gauge

CKD Diagnostic

Preprocessing (Imputation & Encoding)

Model: 100% Recall Classifier

Evaluation & Nephrology Verdict



UNIFIED WEB DEPLOYMENT →  CLINICAL PREDICTION OUTPUT

OUR AI APPROACH

Supervised Learning Architecture

- Algorithmic Paradigm: Supervised learning framework mapping historical clinical parameters directly to validated diagnostic outputs.
- Task Formulation: Formulated as robust classification pipelines predicting discrete medical endpoints and pathological classes.
- Inference Objective: Optimizing high-dimensional decision boundaries to accurately categorize complex physiological signals and patient biometrics.

Clinical Target Taxonomy

- ECG Arrhythmia Telemetry: Multiclass classification mapping 188-sample time-series waveforms across 5 distinct cardiac beat categories.
- Heart Failure Survival Risk: Binary classification evaluating patient survival likelihood versus adverse mortality outcomes.
- Chronic Kidney Disease (CKD): Binary classification distinguishing healthy renal function from pathological CKD presentation across 24 clinical biomarkers.


ECG Heartbeat Classification

Dataset: MIT-BIH Arrhythmia Dataset

Data Type: ECG Signal Data

The dataset contains 5 heartbeat classes

Task: Multiclass Classification



Output:

5 heartbeat classes

Class 0

Class 1

Class 2

Class 3

Class 4

Input:

- ECG heartbeat signal represented by numerical values

ECG – EDA

- **Class Distribution Analysis**
- **Data Visualization**
- **Class Imbalance Detection**

Class 0 → 72,024 (highest class)


Class 1 → 2,223

Class 2 → 5,788

Class 3 → 641 (lowest class)

Class 4 → 6,431

The ECG dataset is highly imbalanced



ECG Model

ECG Preprocessing:

◀ **Feature & Target Separation:**

- $X \rightarrow$ ECG signal features
- $y \rightarrow$ Heartbeat class

◀ **Train-Test Split:**

- 80% Training
- 20% Testing

ECG Data Cleaning:

- Missing values: 0
- Infinite values: 0
- Duplicates: Removed
- Invalid Data

Machine Learning Models

◀ We evaluated:

- Logistic Regression
- K-Nearest Neighbors
- Decision Tree
- Random Forest
- SVM (SVC for Classification)

Model	Main Idea
Logistic Regression	Baseline linear classifier use probability
KNN	Uses nearby samples
Decision Tree	Rule-based decisions
Random Forest	Multiple decision trees
SVM	Finds separating boundaries

	Model	Accuracy	Precision	Recall	F1-Score
0	Logistic Regression	0.673990	0.881088	0.673990	0.737215
1	KNN	0.972730	0.971766	0.972730	0.971559
2	Decision Tree	0.955235	0.954890	0.955235	0.955053
3	Random Forest	0.973004	0.972855	0.973004	0.971147
4	SVM	0.879088	0.897639	0.879088	0.886884

**Best Performing Model:
Random Forest**



Deep Learning

ECG — Deep Learning Models:

- **Neural Network (NN)**
Learns complex relationships between input features
- **Convolutional Neural Network (CNN)**
Extracts important patterns from ECG signals
- **Recurrent Neural Network (RNN)**
Processes sequential ECG data
- **Long Short-Term Memory (LSTM)**
Captures long-term dependencies in ECG sequences




ECG Results

Model Evaluation:

- Recall
- F1-Score
- Confusion Matrix
- Accuracy
- Precision

Best Model:

- Random Forest
 - Accuracy: 97.30%
 - Precision: 97.29%
 - Recall: 97.30%
 - F1-Score: 97.11%



Heart Failure Prediction

Dataset Source & Structure



- Features: 12 Clinical Variables
- Target Variable:
DEATH_EVENT
0 = Survived , 1 = Died

Systemic & Renal Biomarkers




- Serum Creatinine
- Serum Sodium
- Blood Pressure
- Platelets

Heart Function & Enzymes




- Ejection Fraction
- CPK Enzyme

Risk Factors



- Diabetes
- Anemia
- Smoking
- Age Gender
- Follow-up time


Heart EDA



Key Feature Insights:

Time & Ejection Fraction: Strongly linked to higher survival (increases chance of survival)

Serum Creatinine & Age: Highest direct risk factors (increases mortality rate)



Target Distribution

(DEATH_EVENT)

0 (Survived): ~68% (203 patients) — Majority class.

1 (Died): ~32% (96 patients) — Imbalanced target.

Heart Failure Preprocessing



◀ Target Variable: DEATH_EVENT

◀ Outliers Check: Verified extreme values
(Kept as valid medical indicators).

◀ Feature Scaling: Applied StandardScaler to
normalize numeric ranges.

◀ Data Split: Training Set => 80%
Testing Set =>20% (Stratified Split)

Cleaning Data

```
df.isnull().sum()
```

	0
age	0
anaemia	0
creatinine_phosphokinase	0

- Missing Values: 0 missing values
- (Dataset is 100% complete)

```
df.duplicated().sum()  
  
np.int64(0)
```

- Duplicate Records: 0 duplicate rows.



MODELING


=== Machine Learning Models Performance Comparison ===

Model	Accuracy (%)
Gradient Boosting	83.33
Random Forest (Best)	83.33
XGBoost	81.66
Logistic Regression	81.66
AdaBoost	80.00
Bagging	78.33
SVM	76.66
Decision Tree	73.33
KNN	71.66
Naive Bayes	70.00


- **Model Exploration:** Evaluated 10 classification algorithms (Linear, Tree-based, Ensembles).
- **Top Performers:** Gradient Boosting, Random Forest and XGBoost achieved highest accuracies $\geq 81.66\%$.

Model	Key Focus / Characteristics
Logistic Regression	Baseline linear model for direct medical interpretability.
KNN & SVM	Distance & margin boundary classifiers (simpler baseline).
Decision Tree	Single rule-based classifier (prone to overfitting).
Random Forest & Boosting	Ensembles using trees to boost overall accuracy & reduce error.

BEST MODEL



- Random Forest Performance (Best Model)
- Highest Accuracy: Achieved 83.33% overall accuracy.
- Confusion Matrix: Correctly predicted 38 surviving and 12 mortality cases.



- Top Feature Importance
- Primary Driver: time is the most influential variable in prediction.
- Key Indicators: serum_creatinine and ejection_fraction are the 2nd & 3rd top features.

Deep Learning

ANN Model Architecture & Training

- Network Structure: Built a Sequential Artificial Neural Network (ANN) with Dense layers.
- Overfitting Control: Integrated Dropout layers and EarlyStopping to optimize training rounds.
- Optimization: Used Adam optimizer with Binary Crossentropy loss function.
- Test Accuracy: Achieved 85% accuracy on the testing set.

CKD Dataset

- Objective:
 - Early diagnostic classification of Chronic Kidney Disease using multi-modal physiological and biochemical markers.
- Data Hygiene & Schema:
 - Resolved ARFF formatting corruption across 400 patient profiles with 24 clinical features (11 continuous, 13 categorical) and a binary target of 250 CKD and 150 Not-CKD.
- Anti-Leakage Protocol:
 - Executed an 80/20 stratified split (320 train, 80 test) prior to any data transformation, imputation, or statistical analysis.

EDA and Insights

- **Distribution & Class Imbalance**

- Stratified training baseline reflects 200 CKD (62.5%) vs 120 healthy patients (37.5%), preserving exact population proportions before modeling.


- **Biochemical Skew & Waste Accumulation**

- Extreme positive skewness identified in primary metabolic waste indicators: Serum Creatinine (sc = 5.03), Blood Urea (bu = 2.87), and Potassium (pot = 10.51).

- **The Clinical Outlier Insight**

- 100% of IQR outliers across key markers belonged exclusively to CKD patients. Outliers represent biological pathology rather than noise, validating tree-based partitioning.

Target Variable Distribution (CKD vs No CKD)



Model Performance

Key Performance:

- **Top Tier Performance:** Random Forest and AdaBoost tie for the best F1-Score (99.01%), both at 98.75% Accuracy with 100% Recall. No missed CKD cases.
- **Strong Ensembles:** Random Forest, AdaBoost, and SVM all reach 98.75% Accuracy with Random Forest and AdaBoost tying for the top F1-Score (99.01%).
- **Recall Consistency:** Multiple models (Random Forest, Bagging, Gradient Boosting, XGBoost, Decision Tree) achieve 100.0% Recall, ensuring minimal false negatives in disease detection.

Model Selection Rationale

- **Selection Rule:** Since a missed CKD diagnosis is far costlier than a false alarm, we required 100% Recall as a non-negotiable filter first, then used F1-Score as the tiebreaker not accuracy alone, which would be misleading on a 62.5/37.5 class split.
- **Reliable Generalization:** We chose the tuned Random Forest over the tied untuned Random Forest because 5-fold cross-validation gave it a more robust decision boundary (99.51% CV F1), trading a small amount of test F1 (99.01% → 98.04%) for expected stability on unseen patients.
- **Clinical Explainability:** Random Forest preserves transparent feature importance rankings and auditability, allowing clinicians to verify core predictive biomarkers like hemoglobin and serum creatinine.

	Accuracy	Precision	Recall	F1-Score
Random Forest	98.75%	98.04%	100.0%	99.01%
AdaBoost	98.75%	98.04%	100.0%	99.01%
SVM	98.75%	100.0%	98.0%	98.99%
Random Forest (Tuned)	97.5%	96.15%	100.0%	98.04%
Bagging	97.5%	96.15%	100.0%	98.04%
Gradient Boosting	97.5%	96.15%	100.0%	98.04%
XGBoost	97.5%	96.15%	100.0%	98.04%
Logistic Regression	97.5%	98.0%	98.0%	98.0%
Naive Bayes	97.5%	100.0%	96.0%	97.96%
Decision Tree	96.25%	94.34%	100.0%	97.09%
KNN	95.0%	97.92%	94.0%	95.92%

RF Model Evaluation

Holdout Performance (N = 80):

- Overall Test Accuracy: 97.5% (78 / 80 correct classifications).
- True Positives (CKD): 50
- True Negatives (No CKD): 28
- False Positives: 2
- False Negatives: 0

The Zero False-Negative Guarantee:

- Sensitivity / Recall on CKD = 100%: In medical screening, missing a positive patient (Type II error) is catastrophic.
- A 6.7% false alarm rate (2 healthy patients flagged for follow-up testing) is an acceptable clinical trade-off to ensure zero undiagnosed patients slip through

Actual No CKD	28	2
	Predicted No CKD	Predicted CKD
Actual CKD	0	50


Feature Importance

- **Top Predictive Drivers:**

- Hemoglobin (hemo — 0.256): Strongest single split contributor.
- Packed Cell Volume (pcv—0.141): Direct proxy for blood concentration and red cell volume.
- Specific Gravity (sg — 0.127): Urine concentration indicator.
- Serum Creatinine (sc— 0.111): Primary metric of glomerular clearance rate.

- **Biological Grounding:**


- The model learned that impaired kidneys fail to produce erythropoietin, triggering renal anemia (hemo, pcv, rbcc account for ~50% total cumulative importance).
- Structural urine markers (sg, al) combined with waste accumulation (sc) provide orthogonal validation across different stages of kidney damage.



External Validation

- Independent Test Set: The tuned Random Forest was evaluated on a completely separate 200-patient dataset (128 CKD / 72 No CKD) — never seen during training, tuning, or preprocessing fitting.
-
- Performance Held Up: 97.5% Accuracy, 100% Recall, 96.24% Precision, 98.08% F1 which was nearly identical to the original 80-patient holdout test.
- Zero Missed Diagnoses: All 128 CKD cases were correctly identified; the only errors were 5 healthy patients flagged for unnecessary follow-up.
- Takeaway: Consistent performance across two independent test sets indicates the model is learning real clinical signal (hemoglobin, creatinine, specific gravity) rather than memorizing the original 400-patient dataset.

	Metric	Original Test (%)	External Test (%)	Difference (%)
0	Accuracy	97.50	97.50	0.00
1	Precision	96.15	96.24	0.09
2	Recall	100.00	100.00	0.00
3	F1-Score	98.04	98.08	0.05



Deployment Architecture

Framework & Styling



Built with Streamlit utilizing a custom dark-mode UI layout (stApp), radial gradients, and dynamic badge styling for live/demo states.

Modular Navigation



Sidebar routing enables seamless switching between three distinct clinical evaluation engines within a single unified application.

Production Pipeline



Leverages joblib for dynamic loading of pre-trained model binaries (best_ecg_model.pkl, etc.), scalers, and imputers.

Performance Optimization



Utilizes @st.cache_resource for efficient in-memory model caching and rapid asset retrieval during operations.

ECG Model



Input Flexibility

Supports a "Quick Sample Loader" with predefined arrhythmia presets and a "Batch CSV Upload" for raw 188-sample time-series telemetry.



Interactive Telemetry

Renders live Plotly line charts to visualize waveforms, highlighting R-peak detections and mapping raw codes directly into actionable clinical interpretations.


Heart Failure Survival Risk Evaluator

The interface is titled "Manual Patient Assessment" and "Batch CSV Ingestion". It is divided into three main sections: Demographics, Vitals & Lifestyle, and Biomarkers. Each section contains input fields for patient data, including sliders for numerical values and checkboxes for categorical data. A large red button at the bottom is labeled "Evaluate Survival Risk".

Section	Parameter	Value
Demographics	Age	82
	Sex	Female
Vitals & Lifestyle	High Blood Pressure	☐
	Diabetes	☒
	Anaemia	☒
Biomarkers	Ejection Fraction (%)	73
	Serum Creatinine (mg/dL)	7.40
	Serum Sodium (mEq/L)	144
	Creatinine Phosphokinase (mcg/L)	7005
	Platelets (kiloplatelets/mL)	114000

Multi-Tab Interface

Combines "Manual Patient Assessment" (sliders for ejection fraction, creatinine; toggles for comorbidities) with batch CSV ingestion capabilities.



Risk Stratification

Projects continuous mortality probabilities onto a color-coded Plotly Gauge, categorizing patients into Stable, Moderate Caution, or Critical Risk utilizing selective scaling.

CKD Multi-Parameter Diagnostic

- ✓ **Comprehensive Panels** Organizes 24 distinct clinical attributes into logical, collapsible sections: Patient Vitals, Urinalysis, Blood Chemistry, and Clinical History.
- ✓ **Defensive Data Wrangling** Executes robust background preprocessing, handling missing numeric imputation, categorical one-hot encoding, and strict column alignment.
- ✓ **Diagnostic Output** Delivers an immediate clinical verdict (CKD Present vs. Non-CKD / Normal), supported by probability gauges and automated nephrology referral notes.

The screenshot displays two collapsible panels from a diagnostic tool. The 'Urinalysis' panel contains six input fields: Specific Gravity (1.015), Sugar (0), Pus Cell (normal), Bacteria (notpresent), Albumin (0), and RBC (normal). The 'Blood Chemistry' panel contains six input fields: Blood Glucose Random (120), Sodium (138), Packed Cell Volume (40), Blood Urea (40), Potassium (4.40), WBC Count (8000), Serum Creatinine (1.10), Hemoglobin (13.50), and RBC Count (4.80). Each field is represented by a horizontal slider with a red dot indicating the current value.

Panel	Parameter	Value
Urinalysis	Specific Gravity	1.015
	Sugar	0
	Pus Cell	normal
	Bacteria	notpresent
	Albumin	0
	RBC	normal
Blood Chemistry	Blood Glucose Random (mg/dL)	120
	Sodium (mEq/L)	138
	Packed Cell Volume (%)	40
	Blood Urea (mg/dL)	40
	Potassium (mEq/L)	4.40
	WBC Count (cells/cumm)	8000
	Serum Creatinine (mg/dL)	1.10
	Hemoglobin (g/dL)	13.50
	RBC Count (millions/cmm)	4.80

Conclusion and Future Challenges



Core Achievements

Unified Diagnostic Suite: Successfully consolidated complex time-series telemetry and tabular clinical data into a single, low-latency web interface.

Zero-Compromise Sensitivity: Prioritized 100% recall in the CKD classifier to eliminate false negatives, ensuring no critical diagnoses are missed during initial screening.



Future Expansion

Mobile Deployment: Migrating to native mobile apps via Flutter for rapid point-of-care access on physician ward rounds.

Edge Processing: Transitioning real-time ECG waveform classification directly to local microprocessors to reduce server latency.

EHR Integration & Drift: Establishing automated telemetry monitoring and secure ingestion pipelines for live clinical data standards.



Thank you